

# DIACARE 8000 PRO

## TECARTHERAPY

PROFESSIONAL DEVICE TO TAKE CARE OF THE BODY USING THE INTERNAL ENERGIES OF THE BODY AND THE NATURAL AUTOGENERATION MECHANISMS OF TISSUES AND MUSCLES. THE COLOR TOUCH SCREEN NAVIGATION SYSTEM MAKES THE RESEARCH AND THE PARAMETER MODIFICATION EASIER. THE COMPLETE SET OF PROGRAMS, THE ATLAS THAT DESCRIBES WITH PICTURES THE TREATMENTS IN THE DIFFERENT PARTS OF THE BODY AND A SECTION WHERE THE OPERATOR CAN STORE THE TREATMENTS ARE JUST A FEW OF THE REASONS THAT MAKE DIACARE THE BEST EXPRESSION OF GLOBUS TECARTHERAPY.

### TECHNICAL FEATURES

#### Display

Size 7" color touch screen display

#### Power

350W  $\pm$ 10%

#### Outputs

2

#### Emission frequency

400kHz, 448 kHz, 470kHz, 700kHz, 1000kHz, 1200kHz,  $\pm$  10%

#### Treatment modality

Capacitive and resistive

#### Capacitive electrodes

3 size  $\varnothing$  (mm 30, 50, 70)

#### Resistive electrodes

3 size  $\varnothing$  (mm 30, 50, 70)

#### Mains supply

100-240V AC, 50-60Hz, 1.8A

#### Size


mm 515x340x170

#### Weight

22.04 lbs

### EQUIPMENT

- 1 DIACARE 8000 PRO tecartherapy unit
- 1 Large counter electrode and 4 Small (plate)
- 1 Monopolar resistive handpiece with 3 heads  $\varnothing$  (mm 30, 50, 70)
- 1 Monopolar capacitive handpiece with 3 heads  $\varnothing$  (mm 30, 50, 70)
- 1 Bipolar headpiece  $\varnothing$  (mm 50)
- 1 Kit hands free treatment
- 1 Cream can ml 1000
- 1 Power supply unit
- 1 Neutral active return handpiece



CE  
0476



**GLOBUS**  
ITALIAN EXCELLENCE

# DIACARE

## TECARTHERAPY

### PRE-SET PROGRAMS TO AIM THE MOST COMMON PATHOLOGIES

BONE AND ARTICULAR PATHOLOGIES

PATHOLOGIES OF TENDONS AND LIGAMENTS

MUSCULAR AND SOFT TISSUE PATHOLOGIES

VASCULAR AND LYMPHATIC PATHOLOGIES

PATHOLOGIES OF THE NERVOUS SYSTEM



RESISTIVE



CAPACITIVE



BIPOLAR



MULTI-FREQUENCY



PULSED TREATMENTS



TRANSDERMAL VEHICULATION



DISPLAY TOUCH SCREEN



HUMAN ATLAS



HANDS-FREE MODE



PRESET PROGRAMS



AUTOMATIC FREQUENCY SCANNING



ATHERMAL TREATMENTS



COLOR DISPLAY



MAINS SUPPLY